




Via Sant'Anna131-135  
61030 Cartoceto PU (IT)  
Telephone + 39 (0)721830605  
Fax +39 (0)721837154  
e-mail: info@diatheva.com  
www.diatheva.com



## Recombinant Human FKBP12

<b>Description</b>	Recombinant human FKBP12 with N-terminal His-tag, expressed in <i>E. coli</i> and purified through a Nichel-affinity chromatography followed by gel filtration.
<b>Product type</b>	Recombinant protein
<b>Peptide</b>	Full-length protein: Peptidyl-prolyl cis-trans isomerase FKBP1A – Accession number: <b>P62942.2</b>
<b>Expression system</b>	<i>Escherichia coli</i>
<b>Tested by</b>	SDS Page
<b>Purity</b>	>99% as estimated by SDS-PAGE (EU Ph. 5.0 § 2.5.31)

Purified recombinant human FKBP12 protein (lane 1, molecular weight standard; lane 2, 2µg; lane 3, 4µg; lane 4, 6µg) was separated by SDS-PAGE (12% polyacrylamide) and stained with Coomassie Blue.



<b>Form</b>	Liquid
<b>Storage buffer</b>	50mM Hepes pH 8.0, 150mM NaCl, 0,5mM EDTA, 1mM sodium azide.
<b>Storage instructions</b>	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.